



DOWEX UPCORE Mono MA-600

Uniform Particle Size, Macroporous, Strong Base Anion Exchange Resin

| Product | Type | Matrix | Functional group |
|----------------------------|--------------------|--------------------------|------------------|
| DOWEX* UPCORE* Mono MA-600 | Type 1 strong base | Styrene-DVB, macroporous | Quaternary amine |

| Guaranteed Sales Specifications | | Cl ⁻ form |
|-------------------------------------|--|----------------------|
| Total exchange capacity, min. | eq/l kgr/ft ³ as CaCO ₃ | 1.1 24.0 |
| Water content | % | 55 - 65 |
| Bead size distribution [†] | | |
| Mean particle size | µm | 640 ± 50 |
| Uniformity coefficient, max. | | 1.1 |
| >850µ, max. | % | 5 |
| <300µ, max. | % | 0.5 |
| Whole beads, min. | % | 95 |

| Typical Physical and Chemical Properties | | Cl ⁻ form |
|---|----------------------------|----------------------|
| Total swelling (Cl ⁻ → OH ⁻) | % | 15 |
| Particle density | g/ml | 1.06 |
| Shipping weight | g/l lbs/ft ³ | 670 42 |

| Recommended Operating Conditions | |
|----------------------------------|---|
| Maximum operating temperature: | |
| OH ⁻ form | 60°C (140°F) |
| Cl ⁻ form | 100°C (212°F) |
| pH range | 0-14 |
| Bed depth, min. | 1500 mm (4 ft) |
| Pressure drop, design max. | 1.5 bar (22 psi) |
| Pressure drop, max. | 2.5 bar (37 psi) |
| Flow rates: | |
| Service/fast rinse | 5-50 m/h (2-20 gpm/ft ²) |
| Regeneration/displacement rinse | 6-12 m/h (2.4-4.8 gpm/ft ²) |
| Total rinse requirement | 2-4 Bed volumes |
| Regenerant | 2-4% NaOH |
| Organic loading, max. | 5 g KMnO ₄ /l resin |

[†]For additional particle size information, please refer to the Particle Size Distribution Cross Reference Chart (Form No. 177-01775/CH 171-476-E).

*Trademark of The Dow Chemical Company

DOWEX

Ion Exchange Resins

For more information about DOWEX resins,
call Dow Liquid Separations business:

North America 1-800-447-4369
Latin America (+55) 11-5188-9277
Europe (+31) 20-691-6268
Japan (+81) 3-5460-2100
Australia (+61) 2-9776-3226
<http://www.dow.com/liquidseps>

Typical properties and applications:

DOWEX UPCORE Mono MA-600 uniform particle size macroporous strong base anion resin is designed for use in a packed bed counter-current regeneration system. The particle size is specially selected to give a high degree of compaction prior to regeneration and to minimize pressure drop across the bed. The macroporous structure of this resin


makes it the resin of choice in treating waters with high level of organic matter if operating capacity is not of high concern.

Packaging

25 liter or 5 cubic feet fiber drums.

Figure 1. Pressure Drop Data

Temperature = 20° C (68° F)



For other temperatures use:

$$P_T = P_{20^\circ\text{C}} / (0.026 T_{\text{C}} + 0.48), \text{ where } P \equiv \text{bar/m}$$

$$P_T = P_{68^\circ\text{F}} / (0.014 T_{\text{F}} + 0.05), \text{ where } P \equiv \text{psi/ft}$$

Warning: Oxidizing agents such as nitric acid attack organic ion exchange resins under certain conditions. This could lead to anything from slight resin degradation to a violent exothermic reaction (explosion). Before using strong oxidizing agents, consult sources knowledgeable in handling such materials.

Notice: No freedom from any patent owned by Seller or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, Customer is responsible for determining whether products and the information in this document are appropriate for Customer's use and for ensuring that Customer's workplace and disposal practices are in compliance with applicable laws and other governmental enactments. Seller assumes no obligation or liability for the information in this document. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.

Published October 1998.

